

# ArcticZymes Proteinase

ArcticZymes Proteinase is an unspecific endopeptidase originating from an Arctic marine microbial source. It has broad substrate specificity and is easy to inactivate after use.

Histones and other proteins are known to protect nucleic acids from interacting optimally with other DNA binding proteins and enzymes. ArcticZymes Proteinase is ideally suited for transforming chromatin and other dense nucleic acids to naked DNA. The enzyme is easy to heat-inactivate. This allows thermal inactivation at temperatures allowing RNA integrity as well as avoiding dissociation of dsDNA.




Figure 1: ArcticZymes Proteinase is easy to inactivate after use. Figure shows remaining activity after incubation of ArcticZymes Proteinase at different temperatures. Incubation for 30 minutes at indicated temperature in presence of 50 (solid line) or 300 (dotted line) mM NaCl. Remaining activity in % of sample kept on ice in same buffer.

## Unit definition

One unit is defined as the amount of enzyme that produces one  $\mu\text{mol}$  4-nitroaniline (extinction coefficient  $8.8 \text{ mM}^{-1}$  at  $410 \text{ nm}$ ) at  $25^\circ\text{C}$ , pH 8.0 per minute (buffer conditions 50 mM Tris-HCl pH 8.0, 1% DMSO, 1 mM Suc-Ala-Ala-Pro-Phe-pNA).

## Optimal reaction conditions

Optimal enzyme activity is between pH 7–10. Recommended temperature range for activity is from  $25\text{--}40^\circ\text{C}$ . ArcticZymes Proteinase is not dependent of divalent cations such as Calcium ( $\text{Ca}^{2+}$ ) for activity, and is therefore active in buffers containing EDTA ( $<40^\circ\text{C}$  reaction temperature). ArcticZymes Proteinase is compatible with digestion of proteins in the presence of SDS (0.2–1%) and Urea (1–5 M). The enzyme tolerates 500 mM Guanidine Thiocyanate and 1% Triton X-100 ( $>50\%$  activity). No significant loss of activity at 400 mM NaCl is observed.

## Inactivation

Recommended inactivation conditions are incubating the enzyme at 60°C for 15–30 min. Inactivation temperature and time can be adjusted to suit desired application or stringency on residual protease activity. ArcticZymes Proteinase is susceptible to EDTA and DTT at elevated temperatures (>40°C), so EDTA or DTT will give additive effect upon inactivation.

## Inhibitors

The enzyme is inhibited by general serine protease inhibitors, such as PMSE.

## Storage and stability

Optimal stability and storage condition is between pH 6–9. The enzyme is stable upon storage at -20°C for >1 year. ArcticZymes Proteinase is also stable upon storage at 4°C (>6 months) and




Figure 2: The indicated amount of human cells were lysed in 200 µl single cell lysis buffer supplemented with 7.5 U/ml AZ Proteinase (AZProt), using no-proteinase buffer as reference. 1 µl of lysate was used directly in 10 µl qPCR's to quantify human gDNA. AZProt led to improved detection of gDNA for all cell amounts tested.

room temperature (>4 weeks). ArcticZymes Proteinase is stable at dilute concentrations (<0.1 mg/ml) in solutions.

## Ordering information

Product name	Catalogue #	Concentration	Size	Units
ArcticZymes Proteinase	71600-201	200 U/ml	250 µl	50 U
ArcticZymes Proteinase	71600-110	200 U/ml	5 ml	1000 U
ArcticZymes Proteinase	71600-100	According to agreement	According to agreement	According to agreement

## Your OEM partner to deliver novel solutions for genomics and proteomics

ArcticZymes AS  
Sykehusveien 23  
N-9019 Tromsø, Norway

Phone: (47) 7764 8900  
E-mail: [contact@arcticzymes.com](mailto:contact@arcticzymes.com)  
Website: [www.arcticzymes.com](http://www.arcticzymes.com)

ArcticZymes  
987 Old Eagle School Road, Suite 709  
Wayne, PA 19087, USA

Phone: (484) 534 3567  
Fax: (484) 368 3558  
E-mail: [contact-us@arcticzymes.com](mailto:contact-us@arcticzymes.com)

ArcticZymes  
Handelswei 1A 8501 XJ,  
Joure, the Netherlands

Phone: (47) 7764 8900  
E-mail: [contact@arcticzymes.com](mailto:contact@arcticzymes.com)  
Website: [www.arcticzymes.com](http://www.arcticzymes.com)

**Disclaimer:** These products are intended for further manufacturing use or research use only. Certain applications of ArcticZymes AS products may require licenses from others. It is the expressed duty of any receiver of ArcticZymes AS products to acquire such licenses, if necessary. To the extent allowed by law, ArcticZymes AS will not be liable for damages, whether direct, indirect, incidental, or consequential in connection with or arising from this document, including the use of it. ArcticZymes AS products may be covered by pending or issued patents, designs or design applications and/or trademarks or trademark applications or any other registered or unregistered Intellectual Property Right.