

# Shrimp Alkaline Phosphatase

ArcticZymes' Shrimp Alkaline Phosphatase (SAP) is the only heat labile, all-purpose alkaline phosphatase purified from a recombinant source and originally isolated from *Pandalus borealis* (Arctic shrimp). The recombinant production leads to increased storage stability, low batch-to-batch variations and high specific activity. For added flexibility, when lyophilisation may be desired, SAP is also available in a Glycerol FREE format.

## Cloning

SAP offers greater convenience to cloning procedures, since the enzyme may be completely inactivated by a simple heating step. SAP is active in all buffers used for restriction enzymes, and can be added either during restriction digestion, or directly after. With SAP, the user can forget elaborate calculations and multistep incubations, because the enzyme completely dephosphorylates DNA during one, simple, incubation.




Fig 1. Heat inactivation of SAP at 65°C and 75°C

SAP is useful in many molecular biology applications by offering fast and easy dephosphorylation of DNA, RNA and nucleotides. SAP is also active in most restriction enzyme buffers and is completely and irreversibly inactivated after 5 minutes at 65°C. This property simplifies most workflows involving alkaline phosphatase treatment.



### Heat-inactivated

100% heat-inactivated at 65°C



### Enzyme buffers

Works in restriction enzyme buffers



### Clean-up

Used in PCR clean-up



### Removes 5'-phosphates

from DNA, RNA, dNTPs and dephosphorylates proteins



### Glycerol FREE

Also available in Glycerol FREE formulation

## Enzymatic cleanup of PCR products before sequencing





Fig 2. PCR cleanup is important prior to sequencing

Panel A shows a sample that has undergone treatment with SAP and Exonuclease I, while panel B shows a sample without this treatment. Treatment with SAP and Exonuclease I resulted in significant improvement in overall sequence quality.

## Enzymatic cleanup of PCR products before sequencing

Enzymatic cleanup of a PCR product before sequencing eliminates the need for time consuming purification via gels, columns or beads. Simply add SAP and Exonuclease I to your PCR-product and incubate at 37°C for 10 minutes to digest excess primers and nucleotides. Both enzymes are finally inactivated by heating, and the PCR-product is ready for sequencing. The treatment can be done in a single reaction tube, so no further processing is necessary. This enzymatic protocol yields 100% product recovery even for very short PCR products. This PCR cleanup-protocol may also be useful in genotyping and mass spectrometry-based assays.

## Ordering information

	Article no.	Pack Size*	Concentration
<b>SAP</b>	70700-201	1000 U	1 U/μl
	70700-202	5000 U	1 U/μl
	70700-203	25 kU	1 U/μl
	70700-101	200 kU	1 U/μl
	70700-100	Custom	Custom
<b>SAP</b> Glycerol FREE	70710-201	5000 U	≥ 20 U/μl
	70710-100	Custom	Custom

\*One Unit will convert 1 μmol of p-nitrophenyl phosphate per minute to nitrophenol and phosphate at 37°C and pH 10.4 in 0.1 M glycine buffer, 1 mM each of ZnCl<sub>2</sub>, MgCl<sub>2</sub>, and 6 mM 4-nitrophenyl phosphate.

# Your OEM Partner to deliver novel solutions for genomics and proteomics.

## Quality

ArcticZymes is dedicated to the quality of its products and is certified according to ISO 13485:2016. ArcticZymes offers the convenience of providing standard bulk enzymes as off the shelf products. In addition, ArcticZymes offers enzymes in customized formats. For additional information, please contact us.

## Additional information

We are pleased to provide data and information relating to Shrimp Alkaline Phosphatase. Available data includes stability, buffer storage conditions, pH, and specific activity data. For more information, please check our website [www.arcticzymes.com](http://www.arcticzymes.com).

### ArcticZymes Technologies ASA

Sykehusveien 23  
N-9294 Tromsø, Norway

**T** (47) 7764 8900  
**E** [contact@arcticzymes.com](mailto:contact@arcticzymes.com)  
**I** [www.arcticzymes.com](http://www.arcticzymes.com)

### ArcticZymes Technologies ASA

987 Old Eagle School Road, Suite 709  
Wayne, PA 19087, USA

**T** (484) 534 3567  
**F** (484) 368 3558  
**E** [contact-us@arcticzymes.com](mailto:contact-us@arcticzymes.com)

### ArcticZymes AS

Handelswei 1A 8501 XJ  
Joure, the Netherlands

**T** (47) 7764 8900  
**E** [contact@arcticzymes.com](mailto:contact@arcticzymes.com)  
**I** [www.arcticzymes.com](http://www.arcticzymes.com)

These products are intended for further manufacturing use or research use only. Certain applications of ArcticZymes Technologies ASA products may require licenses from others. It is the expressed duty of any receiver of ArcticZymes Technologies ASA products to acquire such licenses, if necessary. To the extent allowed by law, ArcticZymes Technologies ASA will not be liable for damages, whether direct, indirect, incidental, or consequential in connection with or arising from this document, including the use of it. ArcticZymes Technologies ASA products may be covered by pending or issued patents, designs or design applications and/or trademarks or trademark applications or any other registered or unregistered Intellectual Property Right. Version 2.0 • June 2021